SOPIRA® Ject et SOPIRA® Ject vario

- DE Gebrauchsanweisung
GB Instructions for use
FR Mode d'emploi
IT Istruzioni per l'uso

99001540/01

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SOPIRA® Ject und SOPIRA® Ject vario Instructions for use (GB)

Please note:
Due to US Medical Device Legislation and specific requirements for devices marketed in the USA the wording of the user instructions for the US may slightly differ from the standard English version.

Special observance
Please read the following information and application notes carefully. The product may only be used by experienced medical staff in accordance with these instructions. Kulzer does not give any recommendation regarding a method of treatment. The treating medical specialist is responsible for the manner of application and for the selection of the patient. The non-observance or contravention of the instructions will lead to the expiry of the guarantee and to endangering the safety of the patient. If used in combination with further products, the instructions for use and the compatibility statements of those products shall always be taken into account. Please check the product for completeness and intactness before application. If substantiated cause for doubts regarding the completeness and the intactness are on hand, then do not use the product.

Product specification/compatibility
The SOPIRA Ject is available in two various designs. The SOPIRA Ject with swing-out mechanism to open the cartridge holder (the cartridge holder is firmly connected to the syringe), and the SOPIRA Ject vario with slip-on device for separate cartridge holders (the cartridge holders can be replaced). Both models offer self aspiration as well as manual aspiration by means of hooklets.


Indications
Infiltration- and nerve block anaesthesia in dentistry.
Contraindications
The contraindications to be observed are those of the applied anaesthetic. The use of this product is contraindicated in case of known or suspected allergies against components of this product.

Side effects
This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Complications
No complications are known to occur under professional, proper application.

Warning indications
The aspiration hooklets (C) have to be in retracted position when inserting or removing the cartridges. They may otherwise be bent. Extend hooklets only to connect cartridge/plunger, please retract again when separating the connection! Retract the hooklets by turning the knurled screw (H) to the left-hand side until the limiting stop becomes effective!

Special warning indications for use with other compatible products:
– It is essential to clean and sterilize non-sterile medical devices prior to first use!
– Treat a contaminated instrument immediately after each use (see pre-cleaning manual)!
– The instrument must be visually and functionally checked before each use.
– Sort out damaged and defective instruments and replace them.
Only to be used by dentists and for its intended use.

- A Screw thread for cannula.
B Cartridge holder with viewing windows.
C Three pivotable (swing-in and -out) aspiration hooklets at the end of the plunger provide for stable anchoring in the hollow stopper of the cartridge, and they ensure for smooth guidance when injecting and aspirating.
D Joint with swing-out mechanism (for opening the cartridge holder).
E Slip-on device for mounting the cartridge holder.
F Finger rest.
G Plunger.
H Knurled screw to extend and retract the hooklets. The hooklets are extended by turning the knurled screw (H) to the right-hand side (clockwise), and turning left (counter clockwise) will retract them. The screw – mechanism is equipped with a locking device. It is engaged as soon as the hooklets have been completely extended by turning to the right-hand side. An unintentional separation of the connection is impossible. The effect of the locking device may be noticed as a slight initial resistance when turning the knurled screw to the left-hand side to separate the connection.
J Thumb plate fitting.

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Sequence of use
1. Insert the cartridge a) SOPIRA Ject:
Pull the plunger (G) back completely. A spring resistance has to be overcome towards the end. Swing cartridge holder (B) out at a right angle. The holder is in open position. Insert the cartridge. Swing the holder back in alignment with the longitudinal axis. Thereby, completely retract the plunger backwards against the resistance of the spring. Release the plunger against the spring resistance. The cartridge will slip into the correct functioning position, centred in the cartridge holder by a guide ring.

2. Pull the plunger (G) back completely. A spring resistance has to be overcome towards the end. Remove the cartridge holder (B) from the syringe by pulling it out of its rest laterally. Insert the cartridge holder into the cartridge holder. The cartridge holder will slip into the correct functioning position, centred in the cartridge holder by a guide ring.

3. The plunger (G) can be advanced forward several millimetres against the inserted cartridge: The front end of the plunger with the hooklets (C) will be fitted tightly into the recessed void of the hollow stopper. The connection is created by turning the knurled screw (H) to the right-hand side (= extension of the hooklets). The locking device to safeguard the connection is engaged as soon as the hooklets have been completely extended.

4. The cannula is tightly screwed onto the cartridge holder.
The aspiration device is essentially based on the creation of a stable connection between the hollow stopper of the cartridge and the plunger of the syringe, which will bear reactive force and force of pressure. This jam-free fastened connection permits the performance of aspiration and injection without problems. General complications relating to local anaesthesia during dental treatment are often the result of an unintentional intravenous or intra-arterial injection. Failures of onset may also occur in case of an inadvertent intravascular injection, since the local anaesthetic is carried away quickly from the location of the application.

5. The removal of the cartridge is performed in reverse order.

Treatment
General information
1. Whenever you are working with contaminated instruments, follow the personal protection guidelines of the trade association and similar organisations. Wear appropriate protective equipment and ensure that you have had the necessary vaccinations.

2. Risk of infection: Incorrect processing of instruments puts patients, users and third parties at risk of infection and can impair the performance of the instrument.

3. Instruments used on patients known to have or suspected of having Creutzfeldt-Jakob disease or other prion diseases must be disposed of after one use in accordance with specific national requirements.

4. Always comply with the procedures, equipment and devices valid for the user / operator / central sterilisation unit and check them for compatibility with the information provided here.

5. When making up and using solutions, comply with the information on concentration and exposure time provided by the manufacturers of the chemicals. Non-compliance can damage the instrument.

6. Further information on instrument processing can be found at www.k-i.org

Preparation at the place of use
If an instrument has been contaminated, always clean it immediately after use.
To prevent material from drying and adhering to the instrument, large particles of contamination, corrosive solutions and medicinal products must be removed immediately after application of the medicinal product, for example by wiping and rinsing (dry disposal).

Transport
Use suitable transport containers to transport instruments to the reprocessing site in order to rule out the possibility of them constituting a hazard or undergoing external contamination.
Dry disposal is always preferred, whenever possible. Avoid long storage times.

Preparation before mechanical cleaning
If an instrument has been contaminated, always process it immediately after use. If the instrument is a multi-piece device, dismantle it into its constituent parts.
Correct removal of screwed-on cannula, cartridge removal, for dental cylinder cartridge syringes with slip-on device:
Separate the cartridge holder from the base.

Pre-cleaning the surfaces:
Use a brush (no steel brush) or sponge to remove visible contamination or heavy dirt from the surface of the instrument under cold running water (<40°C; drinking water quality).

Pre-cleaning cavities / lumens:
Use a suitable brush (no steel brush) to clean lumens and cavities of the instrument under cold running water (<40°C; drinking water quality). Rinse gaps, grooves and cavities approx. 10 seconds with a pressurised water pistol.

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