**KULZER**  
MITSUI CHEMICALS GROUP

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**GB** - declare under our own responsibility, that the product pressure vessel, according drawing-number (\*), type as above mentioned and serial-number as shown as pressure vessels marking, to which this declaration relates in interpretation, manufacturing, examination and materials, is in conformity with the following standards and normative documents: - declaration of conformity art.(\*\*) according to directive 2009/105 EEC in connection with EC-Type approval acc. art. 10 with the type-approval-N°. (\*\*\*) -

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---

# Operation Instructions

Pressure tank CL-DB for **Kulzer** vacuum pressure casting machines – **Edition: 08-2017**

## Table of contents

Chapter	page
Technical data	7
General requirements	7
Safety information	8
Intended use	9
Connection and operation	9
Maintenance	9
Recurring tests	9
Service	10
Drawing	10
CE Declaration of conformity	10
Document history	11

<b>Technical data</b>	Content:	20 ltr.
	pressure reducer:	0 – 10 bar incl. manometer Ø 40 mm -10 bar
	max. pressure	10 bar
	<b>rated operating pressure:</b>	<b>5 bar</b>
	pressure safety valve:	10 bar
	drainage for condensed water:	manual drain valve at bottom of tank
	air connections:	2 x nozzle Ø 13mm
	accessories:	2 hose clamps (12-25 mm)
	Dimensions:	650 x 200 x 390 mm (without pressure red. & nozzles)

Note: Type designation, series, serial no. and CE – mark can be found on the type plate (at the side of the unit).

## General requirements

Pressure tanks (vessels) may only be operated by persons who are familiar with the handling and can be expected to fulfil their tasks assigned in a reliable manner.

Prior to starting work, the employees must be instructed and receive adequate instructions on

- Operating the pressure tanks
- The special hazards involved in handling pressure tanks and
- Measures to be taken in case of accidents and malfunctions.

Instructions must be provided repeatedly within reasonable intervals. Pressure tanks must be operated in accordance with the intended operating mode.

Pressure tanks must be operated in a way to exclude hazards to employees or third persons, to maintain their operational reliability, to avoid hazardous external corrosion and to protect them against impact loading.

Pressure tanks must be operated in a way to exclude the occurrence of dangerous heating up; the distance to radiators must be at least 0.5 m. Protection against sunlight is not required.

As long as pressure tanks are under pressure, screws of pressure-bearing elements or screwed-in valves must not be removed and may only be tightened by trained staff with suitable tools.

Special incidents, defects and damage to pressure tanks and the equipment of such as well as the activation of their safety devices must be immediately reported to the person responsible for the operation of the tanks.

If a pressure tank reveals any defect or damage which endangers employees or third persons, it must be emptied immediately without any risk; if this is not possible, other

---

suitable measures must be taken which exclude any hazard as far as possible, e.g. evacuating the hazardous area, removing the source of hazard (e.g. compressed air supply).

Pressure tanks must be protected against falling over, falling down, driving against and any impact.

In the event of fire, pressure tanks must be removed from the hazardous area. If this is not possible, the compressed air supply must be disconnected and the pressure tank must be cooled continuously (e.g. by spraying water onto them).

In the event of fire, the fire department must be informed about the presence of pressure tanks.

Pressure tanks which were heated locally or exposed to fire heat require clear labelling and need to be checked prior to any further use, e.g. by trained staff or at the manufacturer's.

## **Safety information**

The pressure tank is subject to the simple pressure vessel directive 2014/29/EU. Any modifications, in particular of the safety fixtures (pressure reducer, manometer, safety valve) are strictly prohibited.

Operator is obliged to observe and fulfil all relevant national standards or safety instructions in accordance with this operating instruction, intended use and safety requirements.

Pressure tanks may only be repaired by trained staff in adequately equipped workshops. Type and extent of repair work must be specified previously.

No welding work or thermal treatment may be performed at the walls of the tank (vessel) exposed to pressure.

The operator must ensure that the internal pressure does not exceed the operating pressure "PS" indicated on the label of the tank (vessel). In case of non-compliance, the manufacturer's operating license will become invalid and explosive pressure losses including the risk of injuries for persons and damage to objects may result.

Vibrating stress which may damage the tank (vessel) as well as corrosion of the tank must be avoided by taking suitable measures. When mounting the tank, it must be ensured that the tank is not screwed firmly to the floor without any vibration absorbers.

**Kulzer** GmbH will not accept any liability for damage to the unit resulting from incorrect operation, non-compliance with these working instructions, inexpert maintenance or repair which has not been performed by the **Kulzer** service agents or if no original spare parts or accessory parts have been used during the exchange of parts.

The design of this pressure tank (BWB Burgau, Germany) is in accordance to simple pressure vessel directive 2014/29/EU. Calculation of design is carried out for static operation condition without any extra load; oscillating load is not considered. The tank is rated for endurance strength until a pressure tolerance of 2 bar (= maximum hysteresis of compressed air pressure).

Any further technical information on the pressure vessel directive 87/404 EU can be found in the publications "Technical rules on the pressure vessel directive (TRB 532)" and "Technical rules for gases" (TRG 280).

The operator is obliged to observe all valid and applicable national standards and safety guidelines etc. about installation and operation of simple pressure vessels. He is obliged to initiate and take all appropriate measures etc. to ensure or maintain the operational safety of the device.

<b>Intended use</b>	<p>The pressure tank must only be used for compressed air in dental laboratories / practice laboratories / universities / dental practices within the scope of technical data. Any other use is not permissible. The device may only be operated in frost-free and dry interior rooms. Access to enable maintenance and to drain condensate must be ensured at <b>any time</b>.</p> <p>The pressure tank serves exclusively as compressed air buffer in the customer's compressed air supply to the vacuum pressure casting machines of <b>Kulzer</b>. The tank is used especially if the diameter or air quantity of the customer's compressed air supply is not sufficient to ensure quick filling of the casting chamber of the vacuum pressure casting machines.</p>
<b>Connection and operation</b>	<p>The compressed air coupling (connection unit) and the hoses used for the vacuum pressure casting machines have been matched with the connections of this tank. Any other couplings or hoses must not be used. Suitable hose sockets for the connection to the customer's compressed air supply can be ordered at Hanau using the Order No. 64600951 for hose, Ø 13 mm (inner) and Order No. 64600945 for hose, Ø 9 mm (inner).</p> <p>The tank should be connected using a direct and thick hose connection (e.g. 13 mm, inner). For this purpose, the hose socket at the tank is used.</p> <p>The operating pressure of the tank should not exceed 5 bar to protect the valves of the casting machines and may have to be adjusted if required.</p> <p>Upon delivery and after connecting it to the compressed air supply and the casting machine, the tank is ready for operation. The operating air pressure of 5 bar needs to be preset at pressure reducer. Pull the black knob up and turn it to adjust the operating pressure. Then pull the knob down to lock.</p> <p><b>Note: Don't turn when knob is locked – danger of damage of pressure reducer!</b></p>
<b>Maintenance</b>	<p>Pressure tanks must be carefully and properly maintained and repaired. Maintenance by a skilled operator includes:</p> <ul style="list-style-type: none"> <li>▪ checking on proper condition, especially on possible leakage</li> <li>▪ checking the proper function of the locking devices (pressure reducer, discharge for condensed water). Locking devices that cannot be unlocked manually must be sent to a specialized company for repair.</li> <li>▪ checking the existence and intactness of the safety devices (safety valve).</li> <li>▪ eliminating minor defects, e.g. tightening of hose clamps.</li> <li>▪ removing contaminations at the outside which may lead to corrosion</li> <li>▪ discharge condensed water and drying the tank</li> </ul> <p>The operator will learn how to perform this maintenance work by reading these operating instructions carefully and completely. For additional repair instructions please contact your local dealer or <b>Kulzer</b> GmbH.</p> <p><b>Instructions for discharge of condensed water:</b></p> <p>Open the discharge bolt at least once a month to check whether condensed water has been accumulated in the tank. In case of regular accumulation of condensed water, the tank must be dried (remove hose connections, heat up by warm air flow to 30 – 40°C, duration: 24 h) and the compressor must be checked and measures for drying compressed air must be taken.</p> <p>Note: Excessive condensation may cause corrosion in the tank and malfunctions in the valve technology of the casting machine.</p>
<b>Recurring Tests</b>	<p>According to the 6. Safety regulation for products for simple pressure tanks and §15 of the Industrial Safety Regulation, recurring tests are required (test of the interior and the strength / stability) after a period of 10 years. The tests need to be performed by qualified staff. The operator is responsible for observing the period specified. Upon expiration of</p>


the period, the approval for use of the tank becomes void. Accordingly, **Kulzer GmbH** has limited the service life of the tank to 10 years. Any further operation after expiration of service life is prohibited and – if not observed - in full and own responsibility of operator. **Please contact us if you have questions concerning qualified disposal!**

**Service**

We appreciate your comments, feedback, and suggestions. **Contact in the countries and more information are available by the displayed QR code or at our website [www.kulzer.com](http://www.kulzer.com)**

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**GB** -declare under our own responsibility, that the product pressure vessel, according drawing-number (\*), type as above mentioned and serial-number as shown as pressure vessels marking, to which this declaration relates in interpretation, manufacturing, examination and materials, is in conformity with the following standards and normative documents: - declaration of conformity art.(\*\*) according to directive 2009/105 EEC in connection with EC-Type approval acc. art. 10 with the type-approval-N°. (\*\*\*)»

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---

## Dokument history

First release:	year 2000
Edition 01.08.03:	accessory kit modified (tube, fittings)
Edition 15.12.03:	incl. pressure reducer, water filter, press. gauge, accessory kit modified
Edition 05.03.04:	drawing updated (position safety valve)
Edition 28.07.04:	additional text for pressure reducer operation / settings
Edition 23.07.15:	Layout change; new company name
Edition 09.08.16:	update of regulatory requirements, standards, guidelines, CE-conformity declaration; operating and service lifetime
<b>Edition 11.08.17:</b>	<b>Adjustments with name changing Kulzer</b>

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Wz./Wu./Noma/Nis/Nosmo/Wj/Wjz/Nisj/Vl, BWB Behälter-Werk Burgau GmbH+Co.KG, StEG, StEG-Strasse 8, D-99331 Burgau -