# Hera

## NEO DUROTERM® 7

DE	Gebrauchsanweisung.....	2
GB	Instructions for use.....	3
FR	Mode d'emploi.....	4
ES	Instrucciones de uso.....	5
IT	Istruzioni per l'uso.....	6
SE	Användningsområden.....	7

EN ISO 7490, Typ 1  
Made in Germany

66009081/11

Manufacturer:  
Kulzer GmbH  
Leipziger Straße 2  
63450 Hanau (Germany)

**KULZER**  
MITSUI CHEMICALS GROUP

# NEO DUROTERM 7

Gebrauchs-  
information

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Instructions  
for use

## Range of applications:

Neo Duroterm 7 is a plaster-bonded, christobalite, nongraphite precision investment material for use when casting precious alloys.

**Mixing ratio:** 100 g powder : 26-30 ml distilled water

**Mixing time, by hand:** 45 seconds

**Mixing time, in a vacuum mixer:** 30 seconds

**Working time:** approx. 5 minutes

**Setting time:** approx. 30 minutes

**Setting expansion:** approx. 0.6 %

**Thermal expansion:** approx. 1.0 %

**Total expansion:** approx. 1.6 %

All data refer to a mixing ratio of 100 g powder : 28 ml distilled water and a test temperature of  $23 \pm 2^\circ\text{C}$ .

## Preheating procedure:

Casting ring				
	x 1	x 3	x 6	x 9
<b>Setting time</b>	30 minutes	30 minutes	30 minutes	30 minutes
<b>Commence preheating</b>	30 minutes	30 minutes	45 minutes	60 minutes
<b>Hold-time at 280°C</b>	30 minutes	30 minutes	40 minutes	50 minutes
<b>Hold-time at 700°C</b>	30 minutes	30 minutes	45 minutes	45 minutes
<b>Heat-rate: max. 7°C/minute</b>				

## Processing:

Weigh out the correct amount of powder and sprinkle it into the correct amount of water. Mix gently with a spatula to wet the powder and then mix in vacuum for 45 seconds.

Apply wetting agent to the wax patterns. Allow the wetting agent on the wax patterns to dry thoroughly before investing. Line the casting ring with non-asbestos liner, leaving at least 3 mm between the liner and the top and bottom edges of the ring.

Due to its low viscosity, Neo Duroterm 7 readily flows around all sections of the wax pattern during investing. After casting, allow the casting rings to cool almost to room temperature. The castings are then easily devested in running water. Due to the fine-grained powder, the castings have smooth surfaces.

The total expansion can be adjusted by varying the mixing ratio of the powder and water.

More water reduces the expansion. Less water increases the expansion.

## Presentation

4.2 kg bag

## Please note:

Store in a dry place.

Neo Duroterm 7 contains quartz powder. Do not inhale the dust.

Technical advice, be it verbal, in writing or by means of tests has been given to the best of our knowledge, but is only to be considered non-binding, also concerning possible rights of third-parties, and does not relieve you of your responsibility to ensure that the products supplied by us are suitable for their intended procedures and purposes. The proper use, application and processing of the products are beyond our control and are, therefore, entirely your responsibility. In the event of Kulzer being liable after all, its liability, regardless of the actual damage, will be limited to the actual cost of the goods supplied by us, and used by you.

We do, of course, guarantee that our products are of perfect quality as stipulated in our general terms of sale and delivery.

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