Hera

LEVOTHERM®

| | | |
|----|----------------------------|----|
| DE | Gebrauchsanweisung | 3 |
| GB | Instructions for use | 5 |
| FR | Mode d'emploi | 7 |
| ES | Instrucciones de uso | 9 |
| IT | Istruzioni per l'uso | 11 |

KULZER
MITSUI CHEMICALS GROUP

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LEVOTHERM is a fast-setting, phosphate-bound precision investment material for use with the model casting technique. After only 30 minutes LEVOTHERM is sufficiently hard for the model to be removed safely from the duplicating material. Models treated with cold hardening liquids are scratch-resistant and have an exceptionally smooth surface. LEVOTHERM is used in conjunction with a special mixing fluid. If distilled water is used, expansion is somewhat more limited. The same liquid, i.e. either water or mixing fluid, should be used for making the model and filling the flask.

Mixing ratio:

| | LEVOTHERM powder | LEVOTHERM liquid |
|----------|------------------|------------------|
| 1 model | 150 g | 21 ml |
| 3 models | 450 g | 63 ml |
| 1 flask | 450 g | 67,5 ml |

Mixing time under vacuum:

60 seconds

Handling time:

3 – 4 minutes, temperature of mixing fluid: 18 – 22°C / 64.4 – 71.6°F

Setting time:

approximately 30 minutes

Time taken for model to dry:

60 minutes at 220°C/ 428°F

Hardening the model:

The model should be a light colour before immersion in the cold hardening liquid for 5 – 8 seconds. After immersion, return the model to the drying chamber for 2 minutes.

Filling the flask:

Mix the required amount of LEVOTHERM and pour down the side of the flask closest to the back of the model. In so doing, hold the flask at a slant righting it again slowly while the mixture is being poured. This process is carried out on a vibrator.

Flask drying:

Heat flask from room temperature to 250°C/482°F and dry initially for 90 minutes at this temperature. Then bring flask up to final temperature.

Initial temperature:

for partial dentures: 950°C/1742°F

for full dentures: 1050°C/1922°F

The flask should be kept at the final temperature for 30 minutes.

Recommended rate of heating to initial temperature: 5°C /41°F per minute.

Removal of model

It is advisable to leave the cast to one side to cool. Where necessary, quenching should take place under a water jet as soon as the red-hot glow has disappeared.

Note on storage of mixing fluid:

The fluid must be protected against frost. Lids of bottles should be replaced after use.

⚠ Safety Instructions:

- Do not open the furnace during preheating as this would cause aggressive combustion.
- The preheat furnace should only be filled to half of its maximum capacity.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer).
Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001. Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before disposing.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

We appreciate your comments, feedback, and suggestions.

Contact in the countries and more information are **available by the displayed QR code** or at our website **www.kulzer.com**

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Kulzer GmbH