Hera

HERAVEST® Saphir

DE	Gebrauchsanweisung	3
GB	Instructions for use	6
FR	Mode d'emploi	9
ES	Instrucciones de uso	12
IT	Istruzioni per l'uso	15
NL	Gebruiksaanwijzing	18
SE	Bruksanvisning	21

KULZER
MITSUI CHEMICALS GROUP

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Graphite-free, phosphate-bonded, precision investment material for „speed“ preheat-ing and programmed preheating when casting precious crown and bridgework.

Storage

Keep the powder in a cool, dry place.

Do not store the mixing liquid at less than +5°C as it is sensitive to low temperatures.

Working temperature

approx. 22°C (Room temperature).

Mixing ratio

100 g Powder – 22 ml Liquid. 60 g Powder – 13 ml Liquid.

Mixing ratio for 160 g Sachet:

Powder 160 g, Liquid 35 ml.

For accurate dispensing please use our dispensing syringes.

Liquid concentration (examples)		
Liquid concentration for 1000 ml liquid		
Concentration	Concentrate (ml)	Distilled water (ml)
50 %	500	500
65 %	650	350
75 %	750	250
95 %	950	50
100 %	1000	–

Concentrations						
Concentration %	Sachet 1 x 160 g		Sachet 2 x 160 g (320 g)		Sachet 3 x 160 g (480 g)	
	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)
100	35	–	70	–	105	–
95	33	2	66	4	99	6
90	32	3	64	6	96	9
85	30	5	60	10	90	15
80	28	7	56	14	84	21
75	26	9	52	18	78	27
70	25	10	50	20	75	30
65	23	12	46	24	69	36
60	21	14	42	28	63	42
55	19	16	38	32	57	48
50	17,5	17,5	35	35	52,5	52,5
45	16	19	32	38	48	57
40	14	21	28	42	42	63

Expansion control (%)						
Indication/ Pattern material		alloy	Casting alloys, high gold content Universal- bonding alloys	Casting alloys reduced gold content Universal bonding alloys	Bonding alloys, high gold content	Bonding alloys reduced gold content
2. Crown and Bridgework	- Wax	65–75	65–75	80–90	80–90	
3. Outer copings	- Pattern resin	70–80	70–80	85–95	85–95	
4. Investment BS Liquid 1 is also available for telescopic and conical crowns. Mixing ratio: 100 g Powder - 22 ml Liquid						

Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

Mixing in vacuum

60 seconds

Place in vacuum for 15 seconds but without switching the paddles on. Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

Working time

Approx. 7 min. at approx. 22°C (Room temperature).

The working time depends on the room temperature. Warmth curtails the working time.

Suitable ring systems	
Use only steel casting rings with cellulose ring liners	Casting ring – X1 and X3 steel: We recommend using one dry cellulose ring liner with Vaseline. Casting ring – X6 and X9 steel: We recommend using two dry cellulose ring liners with Vaseline.

Investing

Fill the ring with investment: only use a vibrator if it is required to improve the flow. Avoid high-vibration! This will lead to the formation of bubbles and breakdown of the mixture.

Placing the rings in a preheated furnace

20 minutes after start of mix. Roughen the top of the mould before placing it in the furnace!

Preheating

Type of alloy	Speed preheating					
	Set-up temp. (°C)	Final temp. (°C)	Hold-times (min.) (commences once the final temp. is attained again)			
			X1	X3	X6	X9
Casting alloys, high gold content	700 – 750	max. 750	30	40	45	60
Universal bonding alloys, high gold content	700 – 800	max. 800	30	40	45	60
Casting alloys, reduced gold content	700 – 750	max. 750	30	40	45	60
Universal bonding alloys, reduced gold content	700 – 800	max. 800	30	40	45	60
Bonding alloys, high gold content	850	max. 950**	30	40	45	60*
Bonding alloys reduced gold content	850	max. 950**	30	40	45	60*

The holding time should be extended if several rings are to be heated at the same time.

* It is not advisable to use X9 casting rings.

** The final temperature can be increased to 950°C if there are casting problems with thin copings. Do not place the mould in the furnace at 950°C, but rather at 850°C, then gradually heat up to the required final temperature at a rate of 10°C per min.

Programmed preheating						
Program level	Temperature (°C)	Heat-rate (°C/min.)	Hold-times (min.) Correct size of mould			
			X1	X3	X6	X9
1	approx. 23°C room temp. - 270°C	8	30	40	45	50
2	270°C - 580°C	8	30	30	30	30
3	580°C - final temp. (950°C)	8	40	40	40	40

⚠ Safety Instructions

- Do not open the furnace during the wax burnout phase.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer).
Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001.
Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

The figures given by our research and development department are guidelines and can be individually adjusted.

www.kulzer.com

We appreciate your comments, feedback, and suggestions.

Contact in the countries and more information are **available by the displayed QR code** or at our website **www.kulzer.com**

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Kulzer GmbH