
MOLDAVEST® master run Hera

DE	Gebrauchsanweisung	3
GB	Instructions for use	6
FR	Mode d'emploi	9
ES	Instrucciones de uso	12
IT	Istruzioni per l'uso	15
RU	Инструкция по применению	18

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Stand: 01.2010

Moldavest® master run Instructions for use GB

Graphite-free, phosphate-bonded, precision investment material for “speed” preheating for investing all types of chrome cobalt molybdenum and precious metal denture bases - for use with silicone duplicating materials.

Storage

Keep the powder in a cool, dry place

Do not store the mixing liquid below +5 °C (41 °F) as it is sensitive to coldness.

Working temperature

approx. 22°C (71.6°F) (room temperature)

Mixing ratio

100 g powder – 18 ml liquid

Mixing ratio for 160 g sachet:

160 g powder – 29 ml liquid

For accurate dispensing please use our dispensing syringes.

Liquid concentration (examples)		
Liquid concentration for 1000 ml liquid		
Concentration	Concentrate (ml)	dist. Water (ml)
50 %	500	500
65 %	650	350
70 %	700	300
75 %	750	250

Concentrations						
	Sachet 1 x 160 g Fabricating 1 Model		Sachet 2 x 160 g (320 g) Fabricating 2 Model		Sachet 3 x 160 g (480 g) Fabricating 3 Model Investing 1 Ring	
Concentration %	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)
90	26	3	52	6	78	9
85	25	4	49	9	74	13
80	23	6	46	12	70	17
75	22	7	44	14	66	21
70	20	9	40	18	60	27
65	19	10	38	20	57	30
60	17	12	34	24	51	36
55	16	13	32	26	48	39
50	14,5	14,5	29	29	43,5	43,5

Expansion control (%)			
Fabricating models		Investing	
Cobalt chrome molybdenum	50–75	Cobalt chrome molybdenum and Gold-Casting Alloys	with Dist. water (H ₂ O)
Gold-Casting Alloys	65		

Stirring

Use a spatula to stir the investment material until the powder has been wetted thoroughly.

Mixing in vacuum

60 seconds Place under vacuum for 15 seconds, but without switching on the mixing paddles. Mix for 60 second. Check frequently to ensure that the vacuum/mixing paddles are functioning correctly. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

Working time

approx. 4 minutes at approx. 22°C (71.6°F) (room temperature)

The working time depends on the room temperature. Warmth curtails the working time.

Fabricating Investment model

Vibrator setting: 6.000 Min.⁻¹

Before pouring the investment into the silicone duplicating mould, spray it with Hera SWE 2000 silicone wetting agent, wait 2 - 3 minutes and blow it thoroughly dry. The excellent flow properties facilitate casting non-porous models. Set the vibrator to high frequency and low intensity. As soon as the duplicating mould has been filled completely, switch the vibrator off. Vibrating too intensely causes the material to segregate, impairs the precision of fit and results in rough surfaces (bubbles) as well as sedimentation.

Drying time Investment model

30 minutes. Remove the model from the duplication mould after 30 minutes.

We recommend applying adhesive to the model or heating the working model briefly to approx. 50°C to ensure that the wax pattern adheres properly.

Full flask investing

Full flask investing Vibrator setting: 6.000 Min.⁻¹

Only mix powder with distilled water (H₂O). Excessive vibration causes sedimentation of the investment and increased bubble formation on the pattern.

Suitable ring systems

- a) Kulzer metal denture base system with paper sleeve
- b) Other alternatives can also be used for casting metal denture bases (e.g. plastic or silicone rings)

Placing the rings

Roughen or grind the top and bottom surfaces of the investment mould before placing it in the furnace.

The casting ring must be placed in the preheating furnace, at final temperature, within 15 minutes of the powder being added to the distilled water (H₂O).

Preheating

Speed preheating			
Type of alloy	Set-up Temp. (°C)	Final Temp. (°C)	Hold-times (min) (commences once the final temp. is attained again.)
Co Cr alloys e.g. Heraenium EH	900–1.030 (note preheating temperature of alloy)	900–1.030 (note preheating temperature of alloy)	45-60
Universal bonding alloys e.g. Mainbond EH	700 (note preheating temperature of alloy)	700 (note preheating temperature of alloy)	60

⚠ Safety Instructions

- Do not open the furnace during the wax burnout phase.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer).
Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001. Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when deinvesting.
- Use an extractor with a fine dust filter when sandblasting.

The data given are furnished to the best of our knowledge and are based on internal test results.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials.

We appreciate your comments, feedback, and suggestions.

Contact in the countries and more information are **available by the displayed QR code** or at our website **www.kulzer-dental.com**

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Kulzer GmbH