Hera

HERAVEST® Speed

DE	Gebrauchsanweisung	3
GB	Instructions for use	6
FR	Mode d'emploi	9
ES	Instrucciones de uso	12
IT	Istruzioni per l'uso	15
NL	Gebruiksaanwijzing	18
SE	Bruksanvisning	21

KULZER
MITSUI CHEMICALS GROUP

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Graphite-free, phosphate-bonded precision investment material for fast preheating for crown and bridgework variably preheatable with and without steel ring.

Storage

Keep the powder in a cool, dry place

Do not store the mixing liquid at less than +5°C (41°F) as it is sensitive to low temperatures.

Working temperature

approx. 22°C (71.6°F) (Room temperature)

Mixing ratio

100 g Powder – 25 ml Liquid

60 g Powder – 15 ml Liquid

Mixing ratio for 160 g Sachet:

Powder 160 g, Liquid 40 ml

For accurate dispensing please use our dispensing syringes.

Liquid concentration (examples)		
Liquid concentration for 1000 ml liquid		
Concentration	Concentrate (ml)	Distilled water (ml)
50 %	500	500
65 %	650	350
75 %	750	250
95 %	950	50

Concentrations						
Concentration %	Sachet 1 x 160 g		Sachet 2 x 160 g (320 g)		Sachet 3 x 160 g (480 g)	
	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)	Liquid (ml)	Dist. water (ml)
100	40	–	80	–	120	–
95	38	2	76	4	114	6
90	36	4	72	8	108	12
85	34	6	68	12	102	18
80	32	8	64	16	96	24
75	30	10	60	20	90	30
70	28	12	56	24	84	36
65	26	14	52	28	78	42
60	24	16	48	32	72	48
55	22	18	44	36	66	54
50	20	20	40	40	60	60
45	18	22	36	44	54	66
40	16	24	32	48	48	72

► Expansion control (%)								
Indication/ Pattern material		alloy		Casting alloys, reduced gold content Universal bonding alloys	Casting alloys, high gold content Universal bonding alloys	Bonding alloys, high gold content	Bonding alloys, reduced gold content	NEM- alloys
		Casting alloys, high gold content Universal bonding alloys	Bonding alloys, high gold content					
1. Inlays, Onlays Partial crowns	- Wax	50 +/-5	50 +/-5	-	-	-	-	-
2. Single crowns and small- span bridges	- Wax	65 +/-5	65 +/-5	75 +/-5	75 +/-5	90 +/-5	90 +/-5	90 +/-5
3. Outer copings telescopic crowns	- Pattern resin	85 +/-5	85 +/-5	90 +/-5	90 +/-5	-	-	-
4. Outer copings tapered crowns	- Pattern- resin	45 +/-5	45 +/-5	-	-	-	-	-
The figures given by our research and development department are guidelines and can be individually adjusted. ► To achieve higher levels of expansion we recommend using Investment Special BS Liquid 3.								

Stirring

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

► Mixing in vacuum

60 seconds

Place in vacuum for 15 seconds but without switching the paddles on. Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

Working time

Approx. 5–6 min. at ap-prox. 22°C (Room temp.)

The working time depends on the room temperature. Warmth curtails the working time.

► Suitable ring systems	
Steel ring: X1 / X3 / X6 / X9	Casting ring–X1 and X3 steel: We recommend using one dry cellulose ring liner with Vaseline. Casting ring–X6 and X9 steel: We recommend using two dry cellulose ring liners with Vaseline.
Ringless: X1 / X3 / X6	Use Kulzer silicone rings for investment.

Investing

Fill the ring with investment: only use a vibrator if it is required to improve the flow. Avoid high vibration! This will lead to the formation of bubbles and breakdown of the mixture.

Placing the rings in a preheated furnace

20 minutes after start of mix. Roughen the top of the mould before placing it in the furnace!

▶ Preheating Speed preheating						
Type of alloy	Set-up temp. (°C)	Final temp. (°C)	Hold-times (min.) (commences once the final temp. is attained again)			
			X1	X3	X6	X9*
Casting alloys, high gold content**	700	max. 700	30	30	45	60
Universal bonding alloys, high gold content**	750-800	max. 800	30	30	45	60
Casting alloys, reduced gold content**	700	max. 700	30	30	45	60
Universal bonding alloys, reduced gold content**	750-800	max. 800	30	30	45	60
Bonding alloys, high gold content	850-900	max. 950	30	30	45	60
Bonding alloys, reduced gold content	850-900	max. 950	30	30	45	60
Non-precious alloys	900-950	max. 950	30	30	45	60

*Fit casting ring at 850°C and heat up to the final temperature of 950°C at maximum heating rate.

Transfer technique:

To achieve better casting results and surface qualities with high-gold content alloys (e.g. Bio Maingold SG) and gold-reduced alloys (e.g. Hera SG) with a preheating temperature of 700°C, first preheat the investment compound to 900°C [depending on casting ring size maintain the temperature (see table)] and then transfer to the 700°C hot oven or reduce the temperature from 900°C to 700°C.

Please follow instructions for the respective alloys.

⚠ Safety Instructions

- Do not open the furnace during the wax burnout phase.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer).
Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001.
Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

⚠ Changes

Revision marking November 2007 / Changes in relation to Issue 02/2005.
Expansion control, vacuum stirring time, suitable casting ring systems, preheating process

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Kulzer equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

The figures given by our research and development department are guidelines and can be individually adjusted.

We appreciate your comments, feedback, and suggestions.
Contact in the countries and more information are available by the displayed QR code or at our website www.kulzer-dental.com

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Kulzer GmbH