# i|BOND® Etch 35 Gel

(DE)	Gebrauchsanweisung.....	3
(GB)	Instructions for use .....	7
(FR)	Mode d'emploi .....	10
(ES)	Instrucciones de uso .....	14
(IT)	Istruzioni per l'uso .....	18
(PT)	Instruções de uso .....	23
(NL)	Gebruiksaanwijzing .....	27
(SE)	Bruksanvisning .....	31
(DK)	Brugervejledning .....	35
(NO)	Bruksanvisning .....	39
(FI)	Käyttöohjeet.....	43
(GR)	Οδηγίες χρήσης .....	47
(CZ)	Návod k použití .....	51
(HU)	Használati utasítás .....	55
(LV)	Lietošanas instrukcija.....	59
(LT)	Naudojimo instrukcija .....	63
(PL)	Instrukcja obsługi .....	67

**KULZER**  
MITSUI CHEMICALS GROUP

### **Produktbeschreibung**

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### **Product description**

iBOND Etch 35 Gel is an etching agent based on orthophosphoric acid. It is used to etch enamel and condition dentine in the context of adhesive filling therapy.

### **Composition**

iBOND Etch 35 Gel contains 35 % phosphoric acid.

The product further contains thickening agents, pigments, water.

### **Indications**

Enamel etching and enamel/dentine etching in the context of the total etch technique prior to the adhesive bonding of:

- Composite and compomer restorations
- Laboratory-fabricated, indirect restorations (e.g. inlays, crowns, bridges, veneers)
- Sealants

### **Contraindications**

- Application to dentine close to the pulp (dentin layer less than 1 mm in thickness)
- Do not use in case of known hypersensitivity to iBOND Etch

### **Side effects**

This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

## **Interactions**

Liner and base materials, such as calcium hydroxide, may be dissolved by the acid, depending on their type and composition.

## **Warnings/precautions**

iBOND Etch 35 Gel contains 35 % phosphoric acid and can cause chemical burns. May irritate eyes. Avoid contact with gingiva, skin, eyes or neighbouring teeth. In case of contact with oral mucosa, rinse thoroughly with plenty of water. If swallowed, thoroughly rinse mouth with water and then drink plenty of water. Consult a physician. In case of contact with skin, thoroughly rinse with soap and water. In the event of contact with eyes, rinse thoroughly with water and seek medical advice. Wear suitable protective gloves. Avoid skin contact. Wear protective clothing. Wear safety glasses/face protection. Immediately remove clothing contaminated with iBOND Etch 35 Gel.

## **Application**

Prior to treatment, carefully clean tooth with fluoride-free polishing paste.

The use of a rubber dam is recommended.

1. After preparation in accordance with adhesive filling therapy, rinse the cavity and dry with oil-free air.
2. Starting from the bevelled enamel areas, apply iBOND Etch to the entire cavity surface, including the dentine. Ensure that the enamel is etched for at least 20 sec and no more than 30 sec. When using iBOND Etch 35 Gel pay special attention to complying with the dwell time on the dentine. Reducing the time to 15 sec in these areas is appropriate to avoid overetching. If desired, only apply iBOND Etch to the enamel and let set for a maximum of 30 sec. Etch unprepared enamel, e.g. in case of fissure sealing, for at least 30 sec and no more than 60 sec. In case of fluoridated enamel, etching for no more than 60 sec may also be required.

**Caution:** In case of deep cavities, we recommend covering the dentine near the pulp with a calcium hydroxide preparation and a base (e.g. glass ionomer cement) prior to etching.

3. Use a water sprayer to completely rinse off iBOND Etch for 20 sec. The cavity must be completely free of any remnants.
4. Dry the cavity with oil-free air (do not dry out dentine!). The etched enamel surface should appear chalky white. In case of contamination of the tooth surface (e.g. with saliva), the etching must be repeated (for a maximum of 10 sec dwell time).
5. Perform further treatment with adhesive (e.g. iBOND® Total Etch) and composite (e.g. VENUS®, CHARISMA®, DURAFILL® VS) in accordance with manufacturer directions.

### **Instructions**

iBOND Etch Gel may be applied using a brush or cannula. Both brushes and cannulas are for one-time use only. To avoid cross infection, items intended for single use must be disposed of immediately after initial use. Release the plunger after applying the gel. Carefully reclose after removing product.

### **Special instructions**

Only to be used by dentists and for its intended use. Store out of reach of children.

### **Storage conditions**

Store at temperatures below 25°C (77°F).

Do not use the material after its expiration date.

Please quote batch number and expiry date in all correspondence about the product.

Expiration date  and batch number : See note on syringe/packaging.

Only recycle completely empty packaging.

Remaining etching gel may be discharged with wastewater if highly diluted with water or neutralised.

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Kulzer GmbH
Kulzer, LLC