DE	Gebrauchsanweisung .....	3
GB	Instructions for use .....	10
FR	Mode d'emploi .....	17
ES	Instrucciones de uso .....	24
IT	Istruzioni per l'uso .....	31
PT	Instruções de uso .....	38
NL	Gebruiksaanwijzing .....	45
SE	Bruksanvisning .....	52
DK	Brugervejledning .....	59
NO	Bruksanvisning .....	66
FI	Käyttöohjeet .....	73

### Produktbeschreibung

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Instructions for use 

## Description

iBOND Self Etch is a light-curing self-priming one-component bonding agent for use in combination with adhesive restorations.

Separate conditioning (etching) of the enamel and dentine is not required (however, use of an additional etching gel on the enamel before application of iBOND Self Etch will not have a negative influence on the bond strength).

iBOND Self Etch was developed for bonding composite resin materials (e.g. composite, compomer, Polyglas®) to the hard tooth structure. iBOND Self Etch etches, primes, bonds and desensitises in one step.

For dental use only. Do not use for indications or applications that are not specifically noted in the instructions for use.

Read the instructions for use carefully before use!  
Follow the safety instructions when using iBOND Self Etch!

## **Composition**

iBOND Self Etch is an acetone/water-based formulation of light-activated methylacrylate resins.

## **Indications**

1. Bonding of direct light-cured composite restorations (including Polyglas® and compomers).
2. Bonding of indirect restorations in combination with a light-curing luting cement: porcelain, Polyglas® and composite restorations (inlays, onlays, veneers, crowns).
3. Sealing hypersensitive areas of teeth.

## **Contra-indications**

Do not use iBOND Self Etch directly on exposed pulp tissue. The use of this product is contraindicated in case of known or suspected allergies against (meth)acrylate compounds. In case of suspected allergies, an allergy test is recommended prior to treatment.

## **Side effects**

This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

## **Interactions**

- The combination of iBOND Self Etch with self-curing composites results in a significant reduction of adhesion and is therefore not recommended.
- Materials containing eugenol may inhibit the polymerisation of iBOND Self Etch.

## **Warnings / Precautions**

- The product is highly flammable.
- The product may irritate the eyes. Avoid contact with eyes. Immediately and thoroughly rinse with plenty of water in case of contact with eyes and contact an ophthalmologist.
- The product contains (meth)acrylate compounds that may cause sensitisation in case of skin contact.
- Wear suitable protective gloves.
- Avoid contact with skin. In case of contact with skin, immediately wash with plenty of water and soap.

## **Application**

### **1. Bonding of direct light-cured composite restorations (including Polyglas® and comonomers).**

#### **1.1. Preparation**

- Clean tooth with an oil and fluoride-free polishing paste.
- Prepare cavity for bonding.
- Rinse cavity with water and air dry.
- Isolate the cavity (the use of a rubber dam is strongly recommended).
- If indicated, place sub-base (e.g. glass ionomer cement).

## 1.2. Dispensing

Bottle application:

- Dispense iBOND Self Etch into the well (1 drop for small cavities, 2 drops for larger cavities).
- Replace bottle cap immediately after dispensing.
- Use iBOND Self Etch immediately after dispensing (no later than 3 minutes).

Single dose application:

- Make sure that the single dose is opened immediately prior to application.
- When removing the liquid rotate the applicator tip or brush briefly in the single dose.

## 1.3. Application

Roughen or grind enamel before application of iBOND Self Etch. Etching/bonding is less effective on unground enamel. On sclerotic dentine separate etching with phosphoric acid for 30 sec is recommended.

- Immediately after dispensing apply a copious amount of iBOND Self Etch to the entire cavity surface and margin with the applicator tip or brush. Make sure that the cavity margin is also sufficiently coated with iBOND Self Etch liquid.
- After application agitate the adhesive slightly for 20 sec. Agitation during the dwell time improves demineralisation and diffusion. Protect the bonding layer from contamination (e.g. by blood or saliva).

- Carefully air-dry iBOND Self Etch with a flow of oil-free air (may take 5–10 sec or longer depending on the geometry of the cavity). The objective is to evaporate the solvent and water from the bonding layer without removing the active ingredients from the tooth surface.

A strong air stream at the start of air-drying will dilute the bonding agent and may result in poor adhesion.

- The surface must be visibly glossy, both after application of iBOND Self Etch and after evaporation of the solvent. The complete cavity surface must be fully coated. If the cavity surface does not appear shiny, apply iBOND Self Etch as described above a second time.
- Light-activate iBOND Self Etch for 20 sec with a halogen or LED light-curing unit. The light-curing time assumes the use of a Kulzer Translux® light-curing unit or other dental curing light of comparable intensity (min. 400–500 mW/cm<sup>2</sup>). If plasma light-curing units (with an output greater than 1200 mW/cm<sup>2</sup>) are used, the curing time can be reduced by 8 sec.

Low-intensity light causes poor adhesion. Curing lights should be checked with reliable light testers at regular intervals. The emitting tip should be placed as close as possible to the curing surface.

- Place restorative material according to manufacturer's instructions.

## **2. Bonding of indirect restorations in combination with a light-curing luting cement: porcelain, Polyglas® and composite restorations (inlays, onlays, veneers, crowns).**

### **2.1. Pre-treatment of restoration**

- Pre-treat the bonding surface of the indirect restoration according to the manufacturer's instructions.

### **2.2. Dispensing and application of iBOND Self Etch**

Light-cure the adhesive prior to application of the luting cement.

- Proceed as described in 1.
- Apply and process luting cement according the manufacturer's instructions. The luting cement must be completely light-cured after insertion of the indirect restoration.

## **3. Sealing hypersensitive areas of teeth.**

### **3.1. Cleaning the tooth**

- Clean tooth with an oil and fluoride-free polishing paste.

### **3.2. Dispensing and application of iBOND Self Etch**

- Proceed as described in 1.
- Carefully remove oxygen-inhibited layer with an alcohol-soaked pellet after light-curing.

If the desensitisation effect is not adequate, apply additional coats of iBOND Self Etch as described, light-cure and remove the oxygen-inhibited layer with a damp pellet.

### **Special Instructions**

- Only to be used by dentists and for its intended use.
- Do not use iBOND Self Etch after the expiry date.

### **Storage**

Do not store above room temperature 25°C (77°F) during daily use. Replace and tighten bottle cap, immediately after use. If not stored correctly the product may lose its effectiveness early. If one of the following is noted, it is a sign that the material is losing its bonding efficacy and should no longer be used:

- Nothing comes out of the bottle when squeezed although there is material left inside.
- Transparent solid or gel-like coagulation is detected in the liquid when dispensed.
- Threads are visible in the material during application with the micro brush.

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Kulzer GmbH
Kulzer, LLC