(DE)	Gebrauchsanweisung	3
(GB)	Instructions for use	11
(FR)	Mode d'emploi	19
(ES)	Instrucciones de uso	27
(IT)	Istruzioni per l'uso	35
(PT)	Instruções de uso.....	44
(NL)	Gebruiksaanwijzing	52
(SE)	Bruksanvisning	60
(DK)	Brugervejledning	68
(NO)	Bruksanvisning	76
(FI)	Käyttöohjeet	84

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Read the directions carefully before use and keep them safe throughout the period the product is in use. Pay attention to the safety information when using iBOND Universal.

Product description

iBOND Universal is a light-curing, self-conditioning all-in-one adhesive for use in adhesive, restorative dentistry.

iBOND Universal was developed for the adhesive bonding and repair of direct filling materials and indirect restorations. Use iBOND Ceramic Primer additionally for the adhesive bonding and repair of glass ceramic.

As a universal adhesive, it may be used as follows depending on the indication:

- SELF-ETCH TECHNIQUE for etching, priming, bonding and desensitisation in a single step.
- SELECTIVE ETCH TECHNIQUE: adhesion to the tooth enamel can be maximised by applying selective enamel etching with phosphoric acid.
- ETCH & RINSE TECHNIQUE (total-etch) with an additional etching step with phosphoric acid on enamel and dentine.

Composition

iBOND Universal is an acetone/water-based solution of light-activated methacrylate monomers.

Indications

1. Bonding of direct restorations for all cavity classes (Black) using light-curing, dual-curing or self-curing methacrylate based composites / compomers.
2. Bonding of light-curing, dual-curing or self-curing core build-up materials.
3. Sealing of hypersensitive tooth areas.
4. Sealing of cavities prior to amalgam restorations.
5. Bonding of fissure sealants.
6. Sealing of cavities and core preparations prior to temporary cementation of indirect restorations (according to the immediate dentine sealing technique).
7. Cementation of indirect restorations with light-curing, dual-curing or self-curing adhesive resin cements.
8. Intraoral repair of composite and compomer restorations, porcelain fused to metal, all ceramic as well as metal restorations.

Contraindications

The use of this product is contraindicated in cases of known or suspected allergies to (meth)acrylate compounds. Do not allow iBOND Universal to come into direct contact with exposed pulp tissue.

Side effects

This product or one of its components may, in particular cases, cause hypersensitive reactions. If suspected, information about the ingredients can be obtained from the manufacturer. In the case of a suspected allergy, we recommend taking an allergy test prior to use.

Interactions

- Materials containing eugenol may impair the polymerisation of iBOND Universal.
- iBOND Universal is light curing and should be protected against ambient light.
- Do not treat the surface of zirconium and aluminium oxide ceramic with phosphoric acid.
- The use of disinfectants, astringents and rinsing solutions containing EDTA or hydrogen peroxide is not recommended. Residues of these materials may impair the adhesion and curing reaction of the adhesive.

But GLUMA Desensitizer and chlorhexidine solutions may be used.

Warning indications

- The product contains acetone and is highly flammable.
- Avoid static charging. Keep away from all sources of ignition.
- Use only in sufficiently well-ventilated rooms. Do not inhale vapours.
- The product may irritate eyes.
- Avoid contact with eyes. Immediately and thoroughly rinse with plenty of water in case of contact with eyes and contact an ophthalmologist.
- Avoid contact with skin. In case of contact with skin, immediately wash with plenty of water and soap.
- Do not apply adhesive directly to the exposed pulp.

Precautions

- Wear appropriate protective clothing, eye protection, a mask and gloves. Eye protection is recommended for the patient. The use of a rubber dam is strongly recommended.
- Seal iBOND Universal and iBOND Ceramic Primer bottles carefully after each use.

- Packages marked with the symbol  are intended for single use. Dispose of after using once and do not use on other patients.
- Do not reuse bottles after contamination with bodily fluids in order to prevent cross-contamination.
- Saliva or blood in the field of application during treatment may lead to failure. Use rubber dam (e.g. Ivory) or other suitable means of isolation.

Special instructions

Only to be used by dentists and for its intended use.

Storage advice

Do not expose to direct sunlight. Store in a sufficiently well ventilated place. Do not store above room temperature 25°C (77°F). Incorrect storage may shorten the shelf-life and the product may prematurely lose its effectiveness.

Do not use after the expiry date. Expiry date and batch number: see note on packaging. Please quote the batch number and expiry date in all correspondence about the product.

Expiry date  and batch number : See note on bottle/packaging.

Disposal instructions

Dispose of in accordance with local regulations.

Application steps for indications 1–7

Preparation of the tooth surface (except 5)

- Prepare cavity according to the rules of the adhesive technique.
- Clean with oil-free and fluoride-free paste.
- Clean with water spray and dry with water-free and oil-free air.

- Apply pulp protection (e.g. calcium hydroxide paste) and/or base (e.g. glass ionomer cement) selectively.

iBOND Universal may be used as a self-etching adhesive *in the self-etch technique*. In combination with untreated tooth enamel, e.g. during fissure sealing (5), after pretreating with polishing paste or powder jets and water spray, etch with 35 % phosphoric acid etching gel (e.g. iBOND Etch 35 Gel) for 20–30 seconds, rinse thoroughly, dry and apply adhesive. Separate etching for 30 seconds is also recommended for sclerotic dentine.

Caution: Pay attention to the directions for using the etching gel.

For *selective enamel etching*, apply 35 % phosphoric acid etching gel (e.g. iBOND Etch 35 Gel) to the prepared enamel and allow to work for 20–30 seconds. Then rinse thoroughly with water and dry with an oil-free air flow. Do not overdry.

In the *etch & rinse technique*, apply 35 % phosphoric acid etching gel (e.g. iBOND Etch 35 Gel) to the prepared enamel and dentine and allow to work (20–30 seconds on enamel/15 seconds on dentine). Then rinse thoroughly with water and dry with an oil-free air flow. Do not overdry.

Preparation of iBOND Universal

From the bottle:

- Drop iBOND Universal into the mixing well (1 drop for small cavities, 2 drops for larger cavities) and use immediately (within 3 minutes).

From the single dose:

- Open the single dose just before use.
- When removing the liquid, briefly rotate the applicator tip or brush in the single dose so that it is completely wetted.

Application

- Use the applicator tip or brush to apply iBOND Universal to the entire surface of the cavity in one generous layer without delay. Take care to ensure that the margin of the cavity is also adequately covered with liquid.
- Gently rub in the adhesive for 20 seconds.
- Air dry carefully with a gentle oil-free air flow until the adhesive film no longer moves. Move the air flow from outside to inside with increasing intensity at the same time.
Caution: If the air flow at the start of blowing is too strong, this leads to the bonding agent thinning out and may result in poor adhesion and post-operative discomfort. Do not leave any excess bonding material/puddles, e.g. on the cavity base.
- The surface must be visibly glossy both after applying iBOND Universal and also after evaporation of the solvent. Ensure complete coverage of the entire cavity surface. If the cavity surface does not appear universally glossy, apply iBOND Universal again as described.

Polymerise iBOND Universal for 10 seconds with a conventional curing lamp ($> 500 \text{ mW/cm}^2$, wavelength range 440–480 nm).

Caution: If the light power is too low this leads to insufficient adhesion. Curing lamps should be inspected at regular intervals with calibrated test equipment. The light emission window should be placed as close as possible to the cavity surface during polymerisation.

Caution: Indication 3 (Sealing hypersensitive areas of teeth.)

Carefully remove oxygen inhibited layer with an alcohol soaked pellet after light curing.

Introduce restoration material immediately in accordance with the manufacturer's instructions. Additional light-curing is recommended for dual-curing cements.

Indication 7 Indirect restorations

Preparation of the material surfaces

Composite and metal restorations in addition to zirconium and aluminium oxide ceramic:

Prepare the restoration in accordance with the manufacturer's instructions. If there is no information available, clean thoroughly with water spray and dry with oil-free air flow.

Caution: Do not etch aluminium and zirconium oxide ceramic with phosphoric acid.

Silicate/glass ceramic and lithium disilicates:

Caution: Ceramic materials require different surface treatments; ask the manufacturer or laboratory about the correct procedure.

Caution: Always etch glass ceramic and IPS Empress surfaces with hydrofluoric acid in accordance with the manufacturer's instructions.

Always apply iBOND Ceramic Primer after pre-treating the ceramic surfaces.

- Measure out the desired quantity of drops of iBOND Ceramic Primer into the mixing well.
- Use the applicator tip or brush to apply to the whole surface of the material without delay and allow to evaporate for 20 seconds.
- Then dry briefly with an oil-free air flow.

Prepare and apply iBOND Universal immediately as previously described (in "Application").

Ensure complete coverage of the entire restoration surface. At the same time take care to ensure that no puddles of adhesive remain before light curing as otherwise the restoration will no longer fit.

After curing of iBOND Universal (see above), apply the adhesive resin cement to the restoration in accordance with the manufacturer's instructions, place on the tooth surface prepared with adhesive and cure.

Indication 8 Intraoral repair

Preparation of the material surfaces

Composite, compomer and metal restorations in addition to zirconium and aluminium oxide ceramic:

Roughen the restorations to be repaired. Rinse thoroughly and dry with an oil-free air flow.

Silicate/glass ceramic and lithium disilicates:

Roughen the surface to be repaired with a diamond bur. Remove the ceramic glaze approx. 1 mm beyond the margin of the fracture site. Rinse and dry with oil-free air. Apply iBOND Ceramic Primer as described previously (Indirect restorations).

Use iBOND Universal as described in "Application". After curing of the adhesive, process appropriate restoration material in accordance with the manufacturer's instructions.

Information which deviates from these directions for use is invalid without the authorisation of Kulzer.

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Kulzer GmbH
Kulzer, LLC