- (DE)** Gebrauchsanweisung
- (GB)** Instructions for use
- (FR)** Mode d'emploi
- (ES)** Instrucciones de uso
- (IT)** Istruzioni per l'uso
- (PT)** Instruções de uso
- (SE)** Bruksanvisning
- (PL)** Instrukcja obsługi

**KULZER**  
MITSUI CHEMICALS GROUP

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## (GB) Instructions for use Pala cre-active®

**Product description**  
 Pala cre-active is a light-curing material for the customised characterisation of methyl methacrylate-based prosthetic resins. It is paste-like and available in a range of colours (white, polar, pink, red, maroon, black, gum, gingiva shade 200, gingiva R50, gingiva light pink, gingiva pink, gingiva clear).

**Composition**  
**Pala cre-active colour fluids**  
 Multifunctional methacrylic acid esters. Silicone dioxide (silanated), photo-initiators, stabilisers, inorganic pigments (some fluorescent).  
 Pala cre-active colour fluids are cadmium-free.

**Pala cre-active gingiva**  
 Multifunctional methacrylic acid esters. Silicone dioxide (silanated), bead polymer, photo-initiators, stabilisers, inorganic pigments.  
 Pala cre-active gingiva is cadmium-free.

**Indications**

- Customisation of (methyl methacrylate-based) total and partial prostheses
- Customisation of implant restorations

**Contra-indications**  
 If the patient is allergic to the one or more ingredients of Pala cre-active, the product must not be used.

**Advice**  
 Pala cre-active gingiva can be mixed with Pala cre-active colour fluids to create customised hues. Only for processing in dental laboratories in accordance with intended use.

**Handling**  
 Prior to the application of Pala cre-active, the surface of the prosthesis is roughened with a coarse diamond bur or via sandblasting (50 µm) and cleaned. Apply the Signum® connector with a brush, allow to set for 2-3 minutes and then polymerise for 90 seconds. The ready-to-use Pala cre-active mixture can then be applied directly to the cured prosthetic resin with a brush or modelling instrument. The different consistencies and hues make it possible to give the different parts of the gums smooth or rough surfaces in pale red or dark

red tones. Each layer is polymerised for 90 seconds.

Pala cre-active can be mixed with Pala cre-active gingiva clear to reduce its colour intensity. If you wish to achieve extremely transparent effects, Pala cre-active gingiva clear can also be used on its own. Following the final polymerisation (180 seconds), the dispersion layer is cleaned off the surface, and the surface is finished off with fine carbide burs, rubber polishers and polishers.

**Polymerisation times (see table below)**

**Storage advice**  
 Do not store above 25°C (77°F). Do not use after expiry date. Avoid direct exposure to sun. Always seal containers again after use. Keep away from children. Visually inspect for damage before first use. Damaged products must not be used.

**Hazard information**  
 Irritates the eyes. Irritates the respiratory system. Irritates the skin. Irritates the mucosa. Keep the container tightly sealed and in a cool place. Keep container in a well-ventilated place. Avoid contact with eyes. In case of contact with eyes, rinse immediately and thoroughly with plenty of water and contact an ophthalmologist. Avoid contact with skin. If product comes into contact with the skin, wash thoroughly with soap and water immediately. Use suitable personal safety equipment (see safety data sheet).

**Disposal information**  
 Regarding disposal, please observe national regulations and the safety data sheet. (Available at www.kulzer-dental.com or from your Kulzer contact.)

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® = registered trademark of Kulzer GmbH