- DE** Gebrauchsanweisung
- GB** Instructions for use
- TR** Kullanma talimatı

Paladent® RR

66051758/11

CE 0197

TR Kullanma talimatı Paladent® RR

Ürün tanımı

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GB Instructions for use Paladent® RR

Product Description

Paladent RR is an autopolymerising product for prosthetic restorations intended for all of the following indications. It consists of a powder component in various colours and a liquid component. Basic prosthetic material for casting permanent and removable dental restorations.

Composition

Cold-curing polymer, consisting of powder and liquid.

Main powder component:

polymethylmethacrylat
Main liquid components: methyl methacrylate, dimethacrylate.

Indications

Dental resin for – completing dentures – restoring the function of dentures, e. g. repairs or relining. For fixed and removable prosthetic devices.

Contra-indications

Direct relining with unpolymerised prosthesis material is contraindicated; for use as intended in the dental laboratory only. A proven allergy to one or more components of Paladent RR is to require urgent medical attention

Notes

It is possible that this product or one of its components may in some cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer. Hypersensitivity reactions (allergies), tingling, unpleasant taste and irritation of the oral mucosa have been described in some cases.

Processing instructions

Hard plaster is recommended for preparing the model. The recommended mixing ratio is 10 g powder:7 mL liquid. The fluid plaster model is insulated twice using the Kulzer insulation material Aislar. Broadly roughen the surfaces of the plastic elements that you are bonding and carefully remove any residual wax. No mechanical retention devices are necessary. Place liquid in mixing container, add the required amount of powder and mix vigourously for approx. 30 seconds. The drippy liquid paste can be worked immediately. The processing time at room temperature 23°C (73°F) is approx. 2 min. Polymerisation in the pressure curing pot (e.g. Palamat®elite) takes 10 min at 55°C (131°F), 2 bar.

Notes

MMA/PMMA-based denture materials always contain a small, permissible amount of residual monomer after fabrication. To further reduce the risk of intolerance reactions, the prosthesis (either newly fabricated or repaired) must be kept in tepid water for at least 12 hours before placement.

Repair

Work completed with Paladent 20 can be repaired and completed using Paladent RR.

Storage advice

Do not store above 25°C (77°F). Do not use after expiry date. Avoid direct exposure to sun. When product has been removed, reseal carefully. Store out of reach of children. Visually inspect for damage before first use. Damaged products must not be used.

Hazard/Safety information

POWDER (dibenzoyl peroxide):

May cause an allergic skin reaction. Avoid breathing dust / fume / gas / mist / vapours / spray. Wear protective gloves / protective clothing / eye protection / face protection. Wash contaminated clothing before reuse. If skin irritation or rash occurs: Get medical advice / attention. IF ON SKIN: Wash with plenty of water.


LIQUID (methyl methacrylate, 2-(2H-Benzotriazol-2-yl)-4-methylphenol):

Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat / sparks / open flames / hot surfaces. - No smoking. Do not breathe dust / fume / gas / mist / vapours / spray. Wear protective gloves / protective clothing / eye protection / face protection. If skin irritation or rash occurs: Get medical advice / attention. In case of fire: Use for extinction: CO₂, sand, extinguishing powder. Store in a well-ventilated place. Keep container tightly closed.




Safety data sheets and more information are available at our website **www.kulzer.com**

® = registered trademark of Kulzer GmbH

Dated: 2017-07

Manufacturer:
Kulzer GmbH
Leipzigzer Straße 2
63450 Hanau (Germany)

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