signum® connector

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GB Instructions for use

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Stand: 2017-07

GB Instructions for use

Light-curing primer for bonding light-curing veneering materials to heat and cold curing PMMA acrylics and artificial dentures.

Composition

Methylmethacrylate, polymethyl-methacrylate, Urethanedimethacrylate, photo-initiators

Range of application

Bonding layer for

- individual shade and shape corrections of denture teeth
- individual shade adaptation of denture base resins

Handling

Roughen the surfaces to be character-ised with a diamond or sandblast with Al₂O₃ (minimum 50 µm at 2 bar pressure). Blow off coarse dust with oil-free air. Fine dust particles that remain on the surface can affect the bonding and should be removed carefully with a clean brush with short, stiff bristles. At this stage, the oil and greasefree surface should not be steam-cleaned or touched. Put some Signum connector in a clean dish, apply a thin layer onto the clean surface with a disposable brush, allow to react for 2–3 mins and then polymerise. Do not touch the surface after polymerisation. If the dispersion layer is wiped off or destroyed by grinding, roughen the sur-face once again, and apply a thin layer of Signum connector. Characterise with Signum cre-active colourfluids and polymerise (refer to the appropriate user instructions). Finally, characterisations are covered with veneering material (dentine, enamel, etc.). Characterisations on denture bases are covered with Palapress®, Palapress® vario, PalaXpress® or PalaXpress® ultra. The PALA gingiva component can also be used as an alternative to Palapress, Palapress vario or PalaXpress denture base materials. All fi nal layers are poly-merised (refer to PALA cre-active user instructions), trimmed and polished to a high gloss in the usual manner.

Product	HiLite power	Heraflash	UniXS®
Signum connector	90 s	90 s	90 s
Signum cre-active/transpafuids	90 s	90 s	90 s
Signum dentine, enamel, effect components	180 s*	180 s*	180 s*

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Advice

Use Signum liquid for bonding cold curing denture base resins to veneers. Apply Signum connector in the same way as previously described. If the reac-tion time is less than 2 mins, the bonding result will be unsatisfactory. Use Signum connector only for Kulzer veneering materials.

Storage advice

Do not store material above 25°C (77°F). Avoid direct sunlight. Always close con-tainers tightly after use. Do not use the material after the expiry date.

Hazard/Safety information

Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respira-tory irritation. Keep away from heat/ sparks/open flames/hot surfaces. – No smoking. Use explosion-proof electrical/ventilating/lighting/equipment. Avoid breathing dust/fume/gas/mist/vapours/ spray. IF ON SKIN (or hair): Remove/Take off immediately all contaminated cloth-ing. Rinse skin with water/shower. Store locked up.

Safety data sheets and more information are available at our website www.kulzer-dental.com

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® = registered trademark of Kulzer GmbH
® = marchio registrato della Kulzer GmbH