

cara Mill

features overview



	cara Mill 2.5	cara Mill 3.5	cara Mill 3.5 PRO	cara Mill 3.5L	cara Mill 3.5L PRO
Processing					
Dry processing	✓✓	✓✓	✓✓	✓✓	✓✓
Wet processing	✓✓	✓✓	✓✓	✓✓	✓✓
Blank processing 98.5 mm	✓✓	✓✓	✓✓	✓✓	✓✓
Block processing (CAD/CAM blocks)	✓✓	✓✓	✓✓	✓✓	✓✓
Materials					
Zirconium oxide/Aluminium oxide	✓✓	✓✓	✓✓	✓✓	✓✓
PMMA/Resins/PEEK/Composite/Wax	✓✓	✓✓	✓✓	✓✓	✓✓
Glass ceramics/Hybrid ceramics	✓✓	✓✓	✓✓	✓✓	✓✓
Sintered metal	✓	✓	✓	✓	✓
Chrome Cobalt	✗	✓	✓✓	✓	✓✓
Titanium	✗	✓	✓✓	✓	✓✓
Applications					
Copings/Crowns/Inlays/Onlays/Veneers	✓✓	✓✓	✓✓	✓✓	✓✓
Bridges up to 16 items	✓✓	✓✓	✓✓	✓✓	✓✓
Telescopes	✓	✓✓	✓✓	✓✓	✓✓
Therapeutic splints (bite splints)	✓✓	✓✓	✓✓	✓✓	✓✓
Tooth models	✓✓	✓✓	✓✓	✓✓	✓✓
Model cast	✓✓	✓✓	✓✓	✓✓	✓✓
Abutment on titanium adhesive basis	✓	✓✓	✓✓	✓✓	✓✓
Prefabricated abutments of Ti/CoCr	✓*	✓✓	✓✓	✓✓	✓✓
One-piece abutments (e. g. PEEK)**	✓	✓✓	✓✓	✓✓	✓✓
One-piece abutments (CoCr, Titanium)**	✗	✓	✓✓	✓	✓✓
Bar structures/One-piece abutment bridges on implants**	✗	✓	✓✓	✓	✓✓

* on demand
 ** advanced knowledge and handling necessary

✓✓ ideally suitable
 ✓ suitable
 ✗ not available

Contact in Germany
 Kulzer GmbH
 Leipziger Straße 2
 63450 Hanau, Germany
 cara-service@kulzer-dental.com
 kulzer.com/cara

Giving a hand to oral health.

KULZER
 MITSUI CHEMICALS GROUP